Spiralbohrer für TCT09 Spiralbohrer für TCT10 Allport
 Twist drills for TCT09 Twist drills for TCT10 Allport
 Taladros espirales para TCT09 Taladros espirales para TCT10 Allport

Schraubendreher
 Screw driver
 Atornillador

CE 0044

Fig.	203RF	203RF	TC203	A2001	HMI41F	31053
Shank ¹	205	205	204	204	205	205
Size	009	011	010	013	023	
Length mm	9.0	9.0	12.0	12.0		

¹204=RA, 205=RAL

²Largest working part diameter in 1/10 mm

Schrauben (je 7 St.)
 Screws (7 p. of each size)
 Tornillos (7 pzs. c/u)

Schrauben (je 7 St.)
 Screws (7 p. of each size)
 Tornillos (7 pzs. c/u)

Fig.	TCT09	TCT09	TCT09	TCT09
Size	007	010	013	016
Length mm	8.0	11.0	14.0	17.0
□	0.9	0.9	0.9	0.9
□	1.1	1.1	1.1	1.1

ONLY in
 Screw System TX Professional

Schrauben (je 3 St.)
 Screws (3 p. of each size)
 Tornillos (3 pzs. c/u)

Schrauben (je 3 St.)
 Screws (3 p. of each size)
 Tornillos (3 pzs. c/u)

Fig.	TCT10	TCT10	TCT10	TCT10
Size	007	010	013	016
Length mm	8.0	11.0	14.0	17.0
□	1.0	1.0	1.0	1.0
□	1.3	1.3	1.3	1.3

ONLY in
 Screw System TX Professional

- Minimaldurchmesser Minimal diameter Diámetro mínimo
- Außendurchmesser External diameter Diámetro externo

Optional erhältlich: Schraubendreher
 Optional available: Screwdriver
 Opción disponible: Atornillador

Art.-No.: TL0T1

Art.-No.: BTX00

Art.-No.: BTXPR

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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
The MEISINGER Screw System TX is an all purpose screw system for easy and safe fixation of bone segments e.g. for horizontal or vertical augmentative grafts. The compatible Torx system enables a secure fixation and an exceptionally high power transmission. In addition, the screws have very flat heads and small diameters (core diameter 0.9 mm or 1.0 mm). In addition to the screws that are in the Screw System TX, the Screw System TX Professional contains screws in the length 13.0 mm and 16.0 mm. Both Screw TX Sets are an ideal addition to all Bone Transfer Systems from our Bone Management® line.

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
To optimize healing, the donor and recipient sites are prepared first with small dotted drillings (blood dots). Then the bone cylinder is placed onto the accordingly prepared bone area and should thereby rest flatly upon the receiving bone. The fixation of the bone cylinder can take place with a hemostat.

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
Subsequently, with the thin spiral drill (1.0 mm) a perforation drilling of the bone cylinder and a drilling up to the counter-cortical of the recipient bone is carried out.

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
In the next step, only the bone cylinder is perforated with the larger twist drill (1.3 mm).

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
Then the screw driver leads the bone screw in the suitable length through the drilling in the bone cylinder up to the bone recipient and sensitively fixes it into the receiving bone area. After a healing process of about 3-4 months the screw can be removed.

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- All purpose screw system for simple and safe fixation of bone segments
- Biocompatible screws made from special titanium alloy in two (resp. four) different lengths
- Both Screw TX Sets are an ideal addition to all Bone Transfer Systems from our Bone Management® line

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General instructions:

- All products delivered are unsterile therefore, before initial and each further potential application to be treated (cleaning / disinfection / sterilisation).

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